



AMWC iCAST GHLS erjekalyon IMCAS ISAPS MEDICA USTEDEN IFATS

AESTHETICS ORTHOBIOLOGICS THERAPEUTICS

T-LAB has been developing products in the field of regenerative medicine since 2012. More than 1 Million+ treatments have been performed with our PRP products; We are Türkiye's leading PRP Kit manufacturer. T-LAB is always up to date with the latest innovations and adopts a people-oriented approach with improved patient outcome as the core of its philosophy, To that end it is closely involved in Regenerative medicine through its Foreign and Local partners, reinforcing its status as a Global Firm and exporting to 65+ countries.

In addition, T-LAB has developed, keeps developing & improving the quality of the treatment that its devices offer through its PRP, PRF, Micrograft Regenerative Kit, Non-Enzymatic SVF Kit, Hydrogel Wound Dressing, and Liposomal Ozone based products. T-LAB always aims to go further in the field of Regenerative Medicine.



REGENERATIVE MEDICINE

DIAQUAL
INTERNATIONAL

DIAQUAL LTD
Upper Mezzanine, 57 Golders
Green Road
London
NW11 8EL



www.diaqual.co.uk | info@diaqual.co.uk | +44 7827 08 2779



Legal Warning:
The content in this brochure has been prepared for general information purposes only and is not intended to provide medical advice on our products, medicine, health or diseases. In addition, the information contained herein does not constitute medical advice regarding medical diagnosis or treatment, nor is it intended to replace the guidance of the physician. Please consult your doctor, before starting any treatment for your disease or health condition regarding the products listed here.

BRO-PRP/ES-EN-R00

DIAQUAL

INTERNATIONAL

CHOICE OF THE BEST

t-lab
PRP
KIT

HIGHER
VOLUME

MORE
PLATELETS

-MORE-
EFFECTIVE
RESULT




UP TO 3 TIMES MORE VOLUME,
MORE PLATELETS

EFFECT OF

3 Classical Session in 1 Session

of T-LAB PRP KIT

CLASSICAL
PRP



Up to
3 times
more
Platelet
Dose

~6cc
PRP

~2cc
PRP

T-LAB
PURE PRP



% 80,9

Platelet Recovery Rate

3.87
BILLION

Platelet Dose

% 95,6

Purity

*Impartial DEPA Analysis Results, IMCAS 2019, Jeremy Magalon, PharmD, PhD, Assoc. Prof.

1 Million+ treatments

Based on several R&D studies implemented within T-LAB's laboratories, it was concluded that the Pure PRP protocol yields "A higher quantity/volume Platelet Recovery Rate" against a classical PRP protocol. As a result, it offers the effect of 3 sessions in one.

BEFORE



AFTER



After 5 PRP sessions
Ali Tuncer, MD



After 3 PRP sessions
Mert Aydoğan, MD



After 3 PRP sessions
Mert Aydoğan, MD



After 1 PRP session
Esin Barak, MD



After 2 PRP sessions
Esin Barak, MD